

CELL AND TISSUE ACQUISITION-DISPENSING SYSTEM

# UNIPICK+ GET A HOLD OF YOUR CELL



The instrument that fits any microscope, collects any cells, and dissects any tissue for any downstream application



### FEATURES







Compatibility with a wide range of inverted microscopes



Full digital control







Y

Collection and deposition of up to 48 samples

# APPLICATIONS

#### Collection and clonal expansion of single cells from adherent cultures



Collection of individual SH-SY5Y (A,B) and CHO (C,D) cells



Clonal expansion of single CHO cell



Viability of collected single cells shown as the percentage of survival rates in dilution control experiments.

### Y

Y

#### Single cell collection and microdissection of native and untreated tissues



Representative collection of a single motor neuron from adult mouse spinal cord. Tissue thickness = 20 μm



Representative microdissecion of fresh frozen adult mouse retina. Tissue thickness =  $12 \mu m$ 



Collection of regions of interest (ROIs) from fixed (archived) tissues



Representative collection of various regions of interest from PFA fixed mouse brain tissue. Tissue thickness = 12  $\mu$ m



# ACCESSORIES



**UNIVERSAL STRADDLE** permits pairing **UNIPICK+** with nearly any inverted microscope from single cell culture to more advanced systems for unmatched flexibility in laboratory work



#### **DISPOSABLE-CAPILLARY UNITS (DCUs)** are ready to use and range in diameter

between 10 to 100  $\mu$ m to ensure desired resolution for the collection of single cells from cultures to microdissection of subanatomical regions from tissue sections



#### UNIVERSAL STAGE ADAPTER

allows for collection and deposition of samples in multiple configurations from slide to slide, plate to plate, slide to plate, slide or plate to PCR tubes, 8-tube strips or a 48-well plate

### SPECIFICATIONS



Interface	Digital
Max Pressure	5 PSI
Max Vacuum	10 PSI
Impulse Duration	0.1 to 1.0 sec
Sample compatibility	Adherent cells, suspended cells,
	3D cultures, native, fresh frozen,
	sucrose treated and fixed tissues
Resolution	Single cell for cultures and
	20 $\mu\text{m}$ ID for tissue sections
Viability of Collected Cells	up to 95%*
Available DCU IDs	<b>10 to 100</b> μ <b>m</b>
Travel Step	<b>1.5</b> μ <b>m</b>
Specimen Holder	25 mm slides,
	35 mm ans 50 mm Petri dishes,
	48 PCR tube holder
Deposition Capacity	up to 48 samples
Acquisition Volume**	DCU20 $\mu m$ - 10 nl to 2.5 $\mu l$
	DCU30 $\mu$ m - 35 nl to 3.0 $\mu$ l
	DCU40 $\mu$ m - 70 nl to 5.0 $\mu$ l
Collection Speed	Tissue sections - up to 12 cells/minute
	Cell culures - up to 25 cells/minute
Microscope compatibility	a wide range of inverted microscopes
	(please, consult NDX representative)

\* - estimated with 3T3, CHO and SH-SY5Y cell lines; \*\* - calculated for standard DMEM media. Acquisition volume may depend on the DCU ID and sample viscosity; \*\*\* - estimated collection speed is given as a reference and may be user/application dependent





NeuroInDx, Inc. 20725 S WESTERN AVE STE 100 TORRANCE, CA 90501-1885



SAFETY INFORMATION

Please read the instruction manual carefully in order to safely operate or handle any NeuroInDx product

Copyright (c) 2016 NeuroInDx, Inc. All rights reserved. NeuroInDx, NeuroInDx logo, UnipicK+ are trademarks of NeuroInDx, Inc. in the U.S. and/or other countries. Trademarks belonging to third parties are the properties of their respective owners. The manufacture, use and/or sale of NeuroInDx, Inc. product(s) may be subject to one or more patents or pending patent applications owned by NeuroInDx, Inc. or licensed to NeuroInDx, Inc. from the third parties. Printed in the USA.

For Research Use Only. Not for use in diagnostic procedures.

### CONTACT US

www.neuroindx.com info@neuroindx.com Tel: 1 424 731 7512 Fax: 1 424 731 7514